




## Parts Included




Drop-In Suction Puck



Screws (4)



Lock Wrench



# CD120SP

## Drop-In Suction Puck

### Fabrication Instructions for 3D printed Sockets



Weight limit: 265 lbs.

2-year warranty against manufacturer defects, excessive wear or breakage.

External Prosthetic Components



EC REP

Advena Limited  
Tower Business Centre  
2nd Flr, Tower Street  
Swatar, BKR 4013  
Malta

CD120SP.revA.11122023

EN   Instructions for Use	PT   Instruções de Utilização
DE   Gebrauchsanweisung	PL   Instrukcja użytkowania
FR   Notice d'utilisation	CS   Návod k použití
ES   Instrucciones para el uso	TR   Kullanım Talimatları
IT   Istruzioni per l'uso	RU   Инструкция по использованию
NO   Bruksanvisning	JA   取扱説明書
DA   Brugsanvisning	ZH   中文说明书
SV   Bruksanvisning	KO   사용 설명서
EL   Οδηγίες Χρήσης	
FI   Käyttöohjeet	
NL   Gebruiksaanwijzing	



www.coyote.us/instructions



Manufactured by Coyote®  
419 N. Curtis Rd., Boise, Idaho 83706 USA  
(208) 429-0026 | www.coyote.us



## Basic 3D Printing Instructions

- 1 Prepare final digital medium for lock attachment with your standard modifications.
- 2 Place digital dummy on distal end of the model in accordance with standard procedure. Digital dummy should be flush with end of socket for best fit.
- 3 Create cavity for drop-in 24.4mm inner height and 62mm diameter. We recommend at least 5 mm socket thickness depending on your printer and materials used.  
*(.stl file is available from Coyote, call 208-429-0026)*
- 4 Print your socket as required. Printing instructions are helpful hints on how to work with the lock and connector. Actual printing thickness and materials are the responsibility of the technician and/or practitioner.

## Finish


- 5 Wipe down O-ring on Drop-In Suction Puck with alcohol before inserting.
- 6 Press lock into place and attach pyramid with supplied 6mm x 18mm screws.
- 7 Attaching pyramid to connector will draw lock into place.
- 8 Use 6x18mm screws provided (see Attention #2 and #4) and Loctite® Blue 242 when attaching pyramid. Torque provided connector screws to 10 Nm.



## Additional Information

### ATTENTION

1. Use the 6x18mm screws provided with typical components. In atypical setups, longer screws may be needed. Always use screws class 10.9 or better.
2. Do not lubricate inside of lock, this will attract debris. If you have a noise issue, it is typically due to seating. Call for technical assistance.
3. Lay-up instructions are helpful hints on how to work with the lock and connector. Actual lay-ups are responsibility of the technician and/or practitioner.
4. If using a flexible inner liner, make sure there is a hole to allow suction. Hole is best place directly over valve entrance. Contact Coyote for more information.
5. Typical Coyote® components use the 6x18mm screws. In atypical setups, longer screws may be needed. Always use screws class 10.9 or better. Make sure screw length fully seats into connector base not just post, longer screws may be needed depending on pyramid thickness. (See image 5a)
6. Always use screws provided during lamination to ensure proper depth is created for attachment.
7. A5. Lay-up instructions are helpful hints on how to work with the lock and connector. Actual lay-ups are responsibility of the technician and/or practitioner.
8. Always use Loctite® Blue 242, Torque settings of connector screws are 10Nm.



## Servicing and Replacing Valve



**9** Insert wrench pins into the holes on the valve body.



**10** Unscrew the valve body counter clockwise.




**11** Remove the valve from the valve body.



**12** At this point, you can blow the valve out with air or rinse it with fresh water, be careful not to lose the screen.



**13** The steel filter disk can be pushed out or blown out. It can be replaced with a new one or cleaned and put back in place. Be careful not to lose it.



**14** Check to make sure the filter and screen are in place. Then, press the valve back into the valve body and screw it into the puck housing.

Patient name: \_\_\_\_\_

For tracking purpose, write LOT number (from funnel of lock) here: \_\_\_\_\_