Parts Included



Housing with attached Connector (*Do not remove Four Hole Connector from lock)







Socket Head Cap Screw (2) (Already installed)



6mm x 18mm Screws (4)

NO| Bruksanvisning

Instructions for Use

Notice d'utilisation

Istruzioni per l'uso

Brugsanvisning

Bruksanvisning

Οδηγίες Χρήσης

Käyttöohjeet

NL | Gebruiksaanwijzing

Gebrauchsanweisung

Instrucciones para el uso

EN I

DF

FR

ES

IT

DA |

SV

FI

FI



Springs (3)





Instruções de Utilização

- Instrukcja użytkowania Návod k použití
- Kullanım Talimatları
- Инструкция по использованию RU
- JA 取扱説明書
- ZH | 中文说明书
- KO | 사용 설명서



www.coyote.us/instructions



Manufactured by Coyote 419 N. Curtis Rd., Boise, Idaho 83706 USA (208) 429-0026 | www.coyote.us





Weight limit: 350 lbs.

2-vear warranty against manufacturer defects, excessive wear or breakage.

Patent No. 6334876 External Prosthetic Components

CD122A3D

3D Air-Lock Drop-In

Fabrication Instructions

for 3D printed Sockets



Coyote

Basic 3D Printing Instructions

Prepare final digital for lock atmedium tachment with your standard modifications.

Place anchor dummy on distal end of socket in accordance with standard procedure for mounting anchor. This will help insure proper lock alianment and depth.

Create cavity for drop-in 49 mm inner height and 61.5 mm diameter. We recommend at least 5 mm socket thickness depending on your printer and materials used. (.stl file is available from Coyote, call 208-429-0026)

Print your socket as required. Printing instructions are helpful hints on how to work with the lock and connector. Actual printing thickness and materials are responsibility of the technician and/or practitioner.

Install Lock into Socket

1

Wipe down O-ring on socket with alcohol before inserting.



attach pyramid with supplied 6mm x 18mm screws.

Attaching pyramid to connector will draw lock into place.

Press lock into place and

Typical Coyote® components use 6x18mm screws provided and Loctite® Blue 242 when attaching pyramid. Torque provided connector screws to 10 Nm. (See Caution #2 and #4)



Carefully smooth inside of hole to allow for easy assembly of lock.



Slide lock plate into lock, springs first. It slides easily ONLY one way. Verify orientation first. (See Caution #3)



Place lock pin in lock to hold lock plate.



Add third spring. Slide 8 release button into valve body.



9 Thread valve body into housing.




Hand-tighten valve 10 body with Coyote Lock Wrench CD103WH or 13mm deep well socket.




Use Coyote Alignment Coupler CD106 for alignment during fitting.

Installing Parts into Air-Lock | See instruction video called "Servicing Air-Lock" at www.coyote.us/airlock




Air-Lock Housing

6



Put your Push Button through the hole on the Valve Body and hand tighten it into the housing.




Put two springs into the two


of the Lock Plate and slide it

side by side circular holes


Hand tighten the Valve Body with a Lock Wrench or 13 mm deep well socket. Tight but not too tight.




Press the Lock Plate with a finger forcing the lock plate into its unlocked position.




Press the push button to compress the springs and this will release your pin.



Push the pin into the funnel hole. This will hold the two springs and lock plate in place.



Installing the parts in the lock housing using this method will help ensure the springs don't aet bent.



Put the third spring into the singular circular hole on the lock plate.

Practitioner Instructions

Poor lock pin spacing leads to premature wear. There should be no play between the lock and liner when fully engaged. You may need to add spacers to the pin to ensure this. Check for proper amount of play before putting lock into socket.




Install pin on liner. Engage lock to check for play between lock and liner.




If there is play, loosen 2 pin away from adapter screw and liner.




Reengage lock to check 3 for play. Repeat until lock seats completely. Remove lock



Gap is created between 4 pin and liner.



5 Based on the gap created by loosening pin, install appropriate number of pin spacers on adapter (see Caution #2)



Replace pin on adapter, making sure base fits snugly on pin spacers



After installing pin spacers, re-engage lock to be sure there is no play.



8

Apply Loctite® Blue 242 to threads of lock pin. Pin may need to be tightened with a 7/16" or 11 mm wrench. (See Caution #4, #5, #12)

Documenting Suction

We view suction not as a component or a code, but as a function. Pistoning and milking can be reduced by maintaining a suction socket when using this lock.

- The suction feature of the lock can be demonstrated and documented very simply.
- Have the amputee step into the lock and seat completely.
- Using the lock wrench, remove the valve body, release button, and outer spring from the lock. The amputee is still locked into the socket, but air is now allowed to flow into the bottom of the socket like a traditional pin.
- Walk the patient normally.
- Amputee may experience a difference in how the socket feels immediately, after some ambulation, or after reinstalling the valve body, release button and outer spring. Patient feedback should be documented.

Call for more information on coding of the Air-Lock: (208) 429-0026.

It is the practitioner's responsibility to demonstrate, document, and select appropriate codes for insurance billing.


Solid Pin install and proper seating Instructions

See instruction video called "CD103P8H Installing Brass Pin" at www.coyote.us/airlock

Poor lock pin spacing leads to premature wear. There should be no play between the lock and liner when fully engaged.




Install pin on liner. Engage lock to check for play between lock and



2 If there is play, loosen pin away from adaptor screw and liner.




Reengage lock to check for play. Repeat until lock seats completely. Remove lock.




If a Gap is created between the pin and liner.

You may need to add spacers to the pin to ensure this. Check for proper amount of play before putting lock into socket.



5 Based on the size of the gap created by loosening pin, install appropriate number of pin spacers on threaded end (see Caution #5).



Replace pin on adaptor, 6 making sure base fits snugly on pin spacers.



After installing pin spacers, re-engage lock to be sure there is no play.



Apply Loctite® Blue 242 8 to threads of lock pin. Pin may need to be tightened with a 7/16" or 11 mm wrench. (See Caution #5)

Air-Lock with P8 Pin

Liner	Size	Spacers used	No. of clicks
Alpha Original	М	1	5
Alpha Select	М	0	5
Ossur	26.5	1	6
Alps	26	1	5


(Chart is a guideline, NOT a guarantee of seating. Verify seating.)

Additional Pins - (in 2 pack - sold separately from lock)




















For tracking purpose, write LOT number (from funnel of lock) here: _



ATTENTION

- Do not position lock with release button pointing posterior or anterior. Typically release button is oriented medially.
- Coyote® Typical components use the atypical 6x18mm screws. ln setups, be Always longer screws may needed. use screws class 10.9 or better. Make sure screw length fully seats into connector base not just post, longer screws may be needed depending on pyramid thickness. Torque provided connector screws to 10 Nm.



- 3. Do not lubricate inside of lock, this will attract debris. If you have a noise issue, it is typically due to seating. Call for technical assistance.
- Always provided during use screws lamination to ensure proper depth is created for attachment.
- Never exceed 3 pin spacers.
- Lay-up and printing instructions are helpful hints on how to work with the lock and connector. Actual lay-ups and printing is the responsibility of the technician and/or practitioner.
- Note number of clicks for engagement. There should be at least 2 to 3 clicks engagement prior to any ambulation and more clicks should occur after a few steps. 5 to 6 clicks (depending on liner) are required for full/ proper seating and engagement.

- Liner threads vary. Begin threading pin into liner by hand whenever possible. A wrench will be needed in cases of tight threads.
- Regardless of threading, always use Loctite® Blue 242 on lock pin threads. Follow liner manufacture instructions as they can vary.
- 10. Do not use the CD103P11 long pin with the Drop-In Easy-Off Lock or the Drop-In Air-Lock. With most liners this longer pin will bottom out in the lock.
- 11. If using a flexible inner liner, do not leave plastic over lock housing, this can cause air leakage and other issues. You should laminate directly over housing. Contact Coyote for more information, or visit the video gallery at www.coyote.us.
- 12. If you have a pin you cannot install, even with a wrench, contact Coyote for a replacement.