




Patient Information

Name: _____

Left Weight: _____

Right Activity Level: _____

Bilateral



AK/Transfemoral Workorder Rev. 04132026

PO#: _____ Date Needed: _____

Received: _____ By: _____

Rush Charge Lab use only

Practitioner: _____ Ship to (if different): _____

Company: _____

Address: _____

City, State, Zip: _____ Shipping Method: _____

Email: _____ UPS/Fedex #: _____

Phone: _____

Prep Fill/Modify are typically 1 unit for pediatric, 2 for adult, and 3 for large AKs/hips/chest

FF Fill _____

FM Modify _____

FA Adjust _____

FT Transfer _____

Ext/Flexion Angle _____

Ab/Adduction Angle _____

Neutral

Connector Alignment Standard offset is 0.5"

Posterior

Anterior

Medial

Lateral

AK Type

FAKT1 AK Test Socket

FAKT2 AK Vivak Test Socket

FAKD1 AK Definitive Socket

Plastic

Copoly

MPE

FHDT 1/4" Heavy Duty Upgrade

Lamination

Light

Standard

FHDL Heavy Duty

Finish

Braid

Solid

Fabric

Pigment

Black

Brown

Caucasian

Latino

Other _____

Inner Flexible

FAF1 AK Ofitras Excel Final Thickness 5/32" 3/16"

FAF2 AK Keasy Cone

None

Brim Type

Full

Gano/Proximal Brim

Add. Padding Add padding to blend out bulbous end

Stovepipe

Valve/Vac

CD105 CQL Valve

Other _____

None

CS Customer Supplied

Quality Assurance

___ Smooth edges ___ Seal test ___ Mods

___ Symmetrical trim ___ Comp. Loctited ___ Pull

___ Pads glued/Incl. ___ All parts Incl. ___ Finish

___ Lock working ___ Proj. tote cleared ___ QA

___ Correct Alignment ___ Engraved P.O. _____ Finish Date

Locks

- CD101 Summit Lock (Retrofit)
- CD102 LowPro Summit
- CD103 Airlock
- CD103D Deep Airlock
- CD103L Lanyard Airlock
- CD103M Medium Airlock
- CD103MD Deep Med. Airlock
- CD103S Small Airlock
- CD103SD Deep Small Airlock
- CD116 Proximal Lock
- CD117 Easy-Off Lock
- CD117D Deep Easy-Off Lock
- CD120 Plastic Lanyard Lock
- CD122A Drop-in Airlock
- CD122E Drop-in Easy-off
- C12S Suction Lanyard Standard
- C12SD Suction Lanyard Drop-In
- Other _____
- None

CS Customer Supplied

-
-
-
-
-
-
-
-
-
-
-
-
-
-
-
-
-
-
-

Connectors/Adaptors

- CD103AF Alignable Four
- CD103MDI Multi-direction Offset
- CD103SDI Single Direction Offset
- CD106 Alignment Coupler
- CD108 Integrator
- CD108S Small Integrator
- CD111 One-Shot Connector (def only)
- CD115CF5 5° Lanyard Connector
- CD119SC Test Socket Connector
- Z3PM Three Prong Male
- Z3PF Three Prong Female
- Z3P Three Prong Thread
- Other _____
- None

CS Customer Supplied


















-
-
-
-
-
-
-
-
-
-
-
-
-
-


Hourly Labor FHL _____


Additional Notes


Send cast/parts to: Coyote, 525 N. Steelhead Boise, ID 83704

TRANSFEMORAL SOCKET MEASUREMENTS

	<u>Anatomical</u>	<u>Reduction</u>	<u>Target</u>	<u>Model</u>	<u>Final Model</u>
M/L					
-Lateral A/P *Standard is 1 Centimeter larger than M/L Mark if different _____					
Proximal Circumference					
Mid Circumference					
Distal Circumference					

-Perineum to distal end 

-Perineum to Floor: 

-KC to Floor: 

Foot Manufacturer _____

Side and Size _____

Gel Liner _____

Lock Mechanism _____

Other Components _____

Send cast and parts to -
Coyote, 525 N. Steelhead Boise, ID 83704