



CD207FK

Fitting Kit

Dynamic Strut *made of Resilergy™*

Weight Range depends on flex - Light to Medium Activity Level



CD207XS - Extra Soft Flex Kit
Weight Range 70-150 lbs.

CD207S - Soft Flex Kit
Weight Range 90-200 lbs.

CD207M - Medium Flex Kit
Weight Range 125-225 lbs.

CD207R - Rigid Flex Kit
Weight Range 200-280 lbs.

CD207XR - Extra Rigid Flex Kit
Weight Range 260-350 lbs.

6-month warranty against manufacturer defects,
excessive wear or breakage.

External Prosthetic Components

CD207FK.post.revA.08142023



www.coyote.us/instructions














Manufactured by Coyote®
419 N. Curtis Rd., Boise, Idaho 83706 USA
(208) 429-0026 | www.coyote.us



Advena Limited
Tower Business Centre
2nd Flr, Tower Street
Swatar, BKR 4013
Malta

Fitting Kit Box for the Dynamic Strut AFO - Map & Parts Included

<p>XSoft, Soft 6x18 (bottom) & (top) 6x20 Medium (Bottom) 6x20</p> 	<p>Dynamic Strut Demos for Fitting - CD207XSS, CD207SS, CD207MS, CD207MS</p> 			
<p>Medium (top) 6x22 Rigid & Extra Rigid 6x22 (bottom) & (top) 6x25</p> 	<p>Dynamic Strut Demos for Fitting - CD207RS, CD207XRS</p> 			
<p>Black Screws 6x14 (8) For XS S 6x16 (8) For M, R, XR</p> 	<p>AFO Tooling Piece (4)</p> 	<p>AFO Upper Mounting Plate (1)</p> 	<p>AFO Lower Mounting Plate (1)</p> 	<p>AFO Plate Foam (4)</p> 
<p>Countersunk Washers (4)</p> 	<p>Shim Button Head Screws XSoft, Soft, Medium 6x20, 6x25 (2 each)</p> 	<p>Shim Self Leveling Washers Washer Bottom (4)</p>	<p>Shim Self Leveling Washers Washer Top (4)</p>	<p>Flexion Shim (2)</p>

EN | Instructions for Use
DE | Gebrauchsanweisung
FR | Notice d'utilisation
ES | Instrucciones para el uso
IT | Istruzioni per l'uso
NO | Bruksanvisning
DA | Bruksanvisning

SV | Bruksanvisning
EL | Οδηγίες Χρήσης
FI | Käyttöohjeet
NL | Gebruiksaanwijzing
PT | Instruções de Utilização
PL | Instrukcja użytkowania
CS | Návod k použití

TR | Kullanım Talimatları
RU | Инструкция по использованию
JA | 取扱説明書
ZH | 中文说明书
KO | 사용 설명서



www.coyote.us/instructions-dynamicstrut