

**Coyote** 

# CD503T

## Amputee Tool Kit



[www.coyote.us/instructions](http://www.coyote.us/instructions)

CD502T.revA.10262022





6 month warranty against manufacturer defects,  
excessive wear or breakage.

Manufactured by Coyote®  
419 N. Curtis Rd., Boise, Idaho 83706 USA  
(208) 429-0026 | [www.coyote.us](http://www.coyote.us)



## Parts Included



### Instructions

- 1 Remove top from base. (*Top is held on by magnets*)
- 2 Put O-rings through side ears on base, wrap O-ring around pylon on prosthetic and wrap O-ring around ears on other side.
- 3 Put 4mm and 5mm hex keys in appropriate vertical slots and attach top with magnet to hold them in place.



#### ATTENTION

This product is for emergency adjustments only. If you need to adjust your prosthetic take it to a qualified technician or prosthetist for any needed adjustments. Be sure to use Loctite® Blue 242 to secure any loose screws and torque to manufacturer specifications.

EN | Instructions for Use  
DE | Gebrauchsanweisung  
FR | Notice d'utilisation  
ES | Instrucciones para el uso  
IT | Istruzioni per l'uso  
NO | Bruksanvisning  
DA | Brugsanvisning  
SV | Bruksanvisning  
EL | Οδηγίες Χρήσης  
FI | Käyttöohjeet  
NL | Gebruiksaanwijzing  
PT | Instruções de Utilização  
PL | Instrukcja użytkowania

CS | Návod k použití  
TR | Kullanım Talimatları  
RU | Инструкция по использованию  
JA | 取扱説明書  
ZH | 中文说明书  
KO | 사용 설명서



[www.coyote.us/instructions](http://www.coyote.us/instructions)