419 N. Curtis Boise, ID 83706 **Patient Information** Coyote Name: 208-429-0026 **Dynamic Strut AFO Workorder** Rev. 11072023 Received: _____By: ____ ☐ Left Weight: Lab use only Practitioner: _____ Ship to (if different): _____ ☐ Right Activity Level: ____ □ Bilateral Address: Prep Fill/Modify are typically 1 unit for pediatric, 2 for adult, and 3 for large AKs/hips/chest City, State, Zip:______Shipping Method:_____ FF Fill ☐ **FM** Modify UPS/Fedex #: ☐ Adjust Phone: No Shell Simple Shell **Molded Shell** No Inner Boot Inner Boot Plastic **Modifications** AFO Type **Thickness** Standard modifications are no reductions on AFOs with 1/8" ankle buildup ☐ **FADS** Dynamic Strut ☐ PolyPro ☐ 3/16" ☐ **1/4" FATHD** Strut Selection 10.5" x 1.5" inches ☐ ProComp[®] **FTPC** ☐ 3/16 ☐ 1/4" CD207XS - Extra Soft Flex Strut Inner Boot Plastic Thickness CD207S - Soft Flex Strut ☐ CD207M - Medium Flex Strut ☐ Yes **FAIB** ☐ No ☐ PolyPro ☐ 3/32" ☐ 1/8" ☐ CoPoly **Anterior Shell Boot Trimline** ☐ MPE 3/4 ☐ No Shell (Strap Only) ☐ Total Contact ☐ FASS - Simple Shell (3/32 MPE) Heel Wedge FAHB (for Stability) Transfer Paper ☐ **FAAS** - Molded Shell Lateral Insole Support | External Posting Medial ZPIS Pelite In- sole (per layer) Straps Ankle Center Calf Color □ Neutral to Bench ☐ Distance from Neutral ____ ZS1 Velcro Black 1"- 1.5" ☐ **ZSHD1** Dacron 1"- 1.5" ☐ Beige Hourly Labor FHL ___ ___ ☐ White Attach Straps **Mod Notes** ☐ Yes ☐ No Add. Padding Thickness Pad Type ☐ Yes **ZAA** ☐ No (Black) (Black) (Skin) ☐ Puff ☐ Plast ☐ P.Cell ☐ Aliplast ☐ Ankle ☐ Puff ☐ Plast ☐ P.Cell ☐ Aliplast ☐ Plantar ☐ Puff ☐ Plast ☐ P.Cell ☐ Aliplast ☐ Calf \square Puff \square Plast \square P.Cell \square Aliplast Additional Notes Other: ☐ Puff ☐ Plast ☐ P.Cell ☐ Aliplast Pad Color ☐ Black ☐ Other: ___ **Quality Assurance** Smooth edges _ Joint Alignment Mods Pull Symmetrical trim Comp. Loctited Finish Pads glued/Incl. ___ All parts Incl. Straps On/Incl. ___ QA Proj. tote cleared

Finish Date | Send cast/parts to: Coyote, 419 N. Curtis Rd., Boise, ID 83706

Strap placement

Engraved P.O.