

**Patient Information**

Name: \_\_\_\_\_



525 N. Steelhead Boise, ID 83704  
208-429-0026

Plastic AFO Workorder Rev. 03252026

Received: \_\_\_\_\_ By: \_\_\_\_\_

- Left Weight: \_\_\_\_\_
- Right Activity Level: \_\_\_\_\_
- Bilateral

PO#: \_\_\_\_\_ Date Needed: \_\_\_\_\_

Rush Charge Lab use only

Practitioner: \_\_\_\_\_ Ship to (if different): \_\_\_\_\_

Company: \_\_\_\_\_

Address: \_\_\_\_\_

City, State, Zip: \_\_\_\_\_

Email: \_\_\_\_\_

Phone: \_\_\_\_\_

Shipping Method: \_\_\_\_\_

UPS/Fedex #: \_\_\_\_\_

**Prep** Fill/Modify are typically 1 unit for pediatric, 2 for adult, and 3 for large AKs/hips/chest

- FF Fill
- FM Modify
- Adjust

**AFO Type**

- FAST Solid Ankle
- FALS Leaf Spring
- FAAT Art. Ankle
- FABVP Boise Bi

**Plastic**

- PolyPro
- CoPoly
- FATHD\* - Heavy Duty
- FTPC ProComp 1/4"

**Thickness**

- 3/32"  1/8"
- 5/32"  3/16"
- 1/4"

**Transfer Paper**

ZTP \_\_\_\_\_

**Inner Boot**

- Yes FAIB  No

**Plastic**

- PolyPro
- CoPoly
- MPE

**Boot Trimline**

- 3/4
- Total Contact

**Thickness**

- 3/32"  1/8"
- 5/32"  3/16"

**Joints**

**Medial**

**Lateral**

- ZTJ Tamarack
- Dorsi \_\_\_\_\_
- Free Motion
- ZBD Double Action
- Spring \_\_\_\_\_
- Pin \_\_\_\_\_

- ZTJ Tamarack
- Dorsi \_\_\_\_\_
- Free Motion
- ZBD Double Action
- Spring \_\_\_\_\_
- Pin \_\_\_\_\_

**90° Stop**

**Heel Buttress FAHB**

- Yes  Lateral
- No  Medial

**Insole Support | External**

- ZPIS Pelite In- \_\_\_\_\_ sole (per layer)
- Neutral to Bench
- Distance from Neutral \_\_\_\_\_

**Straps**

**Ankle Center Calf Color**

- ZS1 Velcro 1"- 1.5" \_\_\_\_\_  Black
- ZSHD1 Dacron 1"- 1.5" \_\_\_\_\_  White
- ZSV1 Velstretch 1"- 2" \_\_\_\_\_

**Attach Straps**

- Yes  No

**Add. Padding**

**Thickness**

**Pad Type**



- Yes ZAA  No
  - Ankle \_\_\_\_\_
  - Plantar \_\_\_\_\_
  - Calf \_\_\_\_\_
  - Other: \_\_\_\_\_
- (Black) (Black) (Skin) (White)
- Puff  Plast  P.Cell  Aliplast
  - Puff  Plast  P.Cell  Aliplast
  - Puff  Plast  P.Cell  Aliplast
  - Puff  Plast  P.Cell  Aliplast
  - Puff  Plast  P.Cell  Aliplast

**Quality Assurance**

- \_\_\_ Smooth edges
- \_\_\_ Symmetrical trim
- \_\_\_ Pads glued/Incl.
- \_\_\_ Straps On/Incl.
- \_\_\_ Strap placement
- \_\_\_ Joint Alignment
- \_\_\_ Comp. Loctited
- \_\_\_ All parts Incl.
- \_\_\_ Proj. tote cleared
- \_\_\_ Engraved P.O.
- \_\_\_ Mods
- \_\_\_ Pull
- \_\_\_ Finish
- \_\_\_ QA
- \_\_\_ Finish Date

**Modifications**

Standard modifications are no reductions on AFOs with 1/8" ankle buildup



**Mod Notes**

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

Hourly Labor  FHL \_\_\_\_\_

**Additional Notes**

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

Send cast/parts to: Coyote, 525 N. Steelhead Boise, ID 83704