

**Patient Information**

Name: \_\_\_\_\_



525 N. Steelhead Boise, ID 83704

208-429-0026

**Dynamic Strut AFO Workorder**

Rev. 03252026

Received: \_\_\_\_\_ By: \_\_\_\_\_

- Left      Weight: \_\_\_\_\_
- Right     Activity Level: \_\_\_\_\_
- Bilateral

PO#: \_\_\_\_\_ Date Needed: \_\_\_\_\_

Rush Charge      Lab use only

Practitioner: \_\_\_\_\_

Ship to (if different): \_\_\_\_\_

Company: \_\_\_\_\_

Address: \_\_\_\_\_

City, State, Zip: \_\_\_\_\_

Shipping Method: \_\_\_\_\_


Email: \_\_\_\_\_

UPS/Fedex #: \_\_\_\_\_

Phone: \_\_\_\_\_

**Prep** Fill/Modify are typically 1 unit for pediatric, 2 for adult, and 3 for large AKs/hips/chest

- FF Fill
- FM Modify
- Adjust



**AFO Type**

**Plastic**

**Thickness**

- FADS Dynamic Strut
- PolyPro
- 3/16"    1/4" FATHD
- ProComp® FTPC    3/16    1/4"

**Inner Boot**

**Plastic**

**Thickness**

- Yes FAIB    No
- PolyPro
- 3/32"
- CoPoly
- 1/8"
- MPE

**Boot Trimline**

- 3/4
- Total Contact

**Heel Wedge**

FAHB (for Stability)

**Transfer Paper**

- Lateral
- Medial
- ZTP \_\_\_\_\_

**Straps**

**Ankle**

**Center**

**Calf**

**Color**

- ZS1 Velcro    1" - 1.5" \_\_\_\_\_
- ZSHD1 Dacron   1" - 1.5" \_\_\_\_\_
- Black
- White

**Attach Straps**

- Yes    No

**Add. Padding**

**Thickness**

**Pad Type**

- Yes ZAA    No
- (Black) (Skin) (Skin) (White)
- Ankle      \_\_\_\_\_
- Plantar    \_\_\_\_\_
- Calf        \_\_\_\_\_
- Other: \_\_\_\_\_
- Puff    Plast    P.Cell    Aliplast
- Puff    Plast    P.Cell    Aliplast
- Puff    Plast    P.Cell    Aliplast
- Puff    Plast    P.Cell    Aliplast

**Quality Assurance**

- \_\_\_\_ Smooth edges      \_\_\_\_ Joint Alignment      \_\_\_\_ Mods
- \_\_\_\_ Symmetrical trim    \_\_\_\_ Comp. Loctited        \_\_\_\_ Pull
- \_\_\_\_ Pads glued/Incl.      \_\_\_\_ All parts Incl.        \_\_\_\_ Finish
- \_\_\_\_ Straps On/Incl.        \_\_\_\_ Proj. tote cleared    \_\_\_\_ QA
- \_\_\_\_ Strap placement      \_\_\_\_ Engraved P.O.        \_\_\_\_ Finish Date

**Modifications**

Standard modifications are no reductions on AFOs with 1/8" ankle buildup

**Strut Selection**

Standard 10.5" x 1.5" inches / Mini 8.5"x1.25"

**Standard Dynamic Strut**

- CD207XS - Extra Soft Flex Strut
- CD207S - Soft Flex Strut
- CD207M - Medium Flex Strut
- CD207R - Rigid Flex Strut
- CD207XR - Extra Rigid Flex Strut

**Mini Dynamic Strut**

- CD209S - Soft Flex Strut
- CD209M - Medium Flex Strut
- CD209R - Rigid Flex Strut
- CD209XR - Extra Rigid Flex Strut

**Anterior Shell**

- No Shell (Strap Only)
- FASS - Simple Shell
- FAAS - Molded Shell

**Insole Support / External Posting**

- ZPIS Pelite Insole (per layer)
- Neutral to Bench
- Distance from Neutral \_\_\_\_\_

Hourly Labor  FHL \_\_\_\_\_

**Mod Notes**

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

**Additional Notes**

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

Send cast/parts to: Coyote, 525 N. Steelhead Boise, ID 83704