

Patient Information

Name: _____



419 N. Curtis Boise, ID 83706

208-429-0026

Dynamic Strut AFO Workorder

Rev. 01162024

Received: _____ By: _____

- Left Weight: _____
- Right Activity Level: _____
- Bilateral

PO#: _____ Date Needed: _____

Rush Charge Lab use only

Practitioner: _____

Ship to (if different): _____

Company: _____

Address: _____

City, State, Zip: _____

Shipping Method: _____


Email: _____

UPS/Fedex #: _____

Phone: _____

Prep Fill/Modify are typically 1 unit for pediatric, 2 for adult, and 3 for large AKs/hips/chest

- FF Fill
- FM Modify
- Adjust



AFO Type

Plastic

Thickness

- FADS Dynamic Strut PolyPro 3/16" 1/4" FATHD
- ProComp® FTPC 3/16 1/4"

Inner Boot

Plastic

Thickness

- Yes FAIB No
- PolyPro 3/32" 1/8"
- CoPoly
- MPE

Boot Trimline

- 3/4
- Total Contact

Heel Wedge

FAHB (for Stability)

Transfer Paper

- Lateral
- Medial

ZTP _____

Straps

Ankle

Center

Calf

Color

- ZS1 Velcro 1" - 1.5" _____ _____ _____ Black
- ZSHD1 Dacron 1" - 1.5" _____ _____ _____ Beige
- _____ _____ _____ _____ _____ White

Attach Straps

- Yes No

Add. Padding

Thickness

Pad Type

- Yes ZAA No
- (Black) (Black) (Skin) (White)
- Ankle _____ Puff Plast P.Cell Aliplast
- Plantar _____ Puff Plast P.Cell Aliplast
- Calf _____ Puff Plast P.Cell Aliplast
- Other: _____ _____ Puff Plast P.Cell Aliplast

Pad Color

- Black Other: _____

Quality Assurance

- ____ Smooth edges ____ Joint Alignment ____ Mods
- ____ Symmetrical trim ____ Comp. Loctited ____ Pull
- ____ Pads glued/Incl. ____ All parts Incl. ____ Finish
- ____ Straps On/Incl. ____ Proj. tote cleared ____ QA
- ____ Strap placement ____ Engraved P.O. ____ Finish Date

Modifications

Standard modifications are no reductions on AFOs with 1/8" ankle buildup

Strut Selection

Standard 10.5" x 1.5" inches / Mini 8.5"x1.25"

Standard Dynamic Strut

- CD207XS - Extra Soft Flex Strut
- CD207S - Soft Flex Strut
- CD207M - Medium Flex Strut
- CD207R - Rigid Flex Strut
- CD207XR - Extra Rigid Flex Strut

Mini Dynamic Strut

- CD209S - Soft Flex Strut
- CD209M - Medium Flex Strut
- CD209R - Rigid Flex Strut

Anterior Shell

- No Shell (Strap Only)
- FASS - Simple Shell (3/32 MPE)
- FAAS - Molded Shell

Insole Support / External Posting

- ZPIS Pelite Insole (per layer)
- Neutral to Bench
- Distance from Neutral _____

Hourly Labor FHL _____

Mod Notes

Additional Notes

Send cast/parts to: Coyote, 419 N. Curtis Rd., Boise, ID 83706