

**Patient Information**

Name: \_\_\_\_\_



525 N. Steelhead Boise, ID 83704

208-429-0026

**Dynamic Strut AFO Workorder**

Rev. 10142025

Received: \_\_\_\_\_ By: \_\_\_\_\_

- Left      Weight: \_\_\_\_\_
- Right     Activity Level: \_\_\_\_\_
- Bilateral

PO#: \_\_\_\_\_ Date Needed: \_\_\_\_\_

Rush Charge      Lab use only

Practitioner: \_\_\_\_\_

Ship to (if different): \_\_\_\_\_

Company: \_\_\_\_\_

Address: \_\_\_\_\_

City, State, Zip: \_\_\_\_\_

Shipping Method: \_\_\_\_\_


Email: \_\_\_\_\_

UPS/Fedex #: \_\_\_\_\_

Phone: \_\_\_\_\_

**Prep** Fill/Modify are typically 1 unit for pediatric, 2 for adult, and 3 for large AKs/hips/chest

- FF Fill
- FM Modify
- Adjust



- AFO Type**      **Plastic**      **Thickness**
- FADS Dynamic Strut     PolyPro       3/16"     1/4" FATHD
  - ProComp® FTPC     3/16       1/4"

- Inner Boot**      **Plastic**      **Thickness**
- Yes FAIB     No       3/32"
  - CoPoly       1/8"
  - MPE

- Boot Trimline**
- 3/4
  - Total Contact

- Heel Wedge**    FAHB (for Stability)      **Transfer Paper**
- Lateral
  - Medial
  - ZTP \_\_\_\_\_

- Straps**      **Ankle**    **Center**    **Calf**    **Color**
- ZS1 Velcro    1"- 1.5"    \_\_\_\_\_    \_\_\_\_\_    \_\_\_\_\_     Black
  - ZSHD1 Dacron 1"- 1.5"    \_\_\_\_\_    \_\_\_\_\_    \_\_\_\_\_     Beige
  - \_\_\_\_\_    \_\_\_\_\_    \_\_\_\_\_    \_\_\_\_\_    \_\_\_\_\_     White

- Attach Straps**
- Yes     No

- Add. Padding**    **Thickness**      **Pad Type**
- Yes ZAA     No
  - Ankle      \_\_\_\_\_       Puff  Plast  P.Cell  Aliplast
  - Plantar    \_\_\_\_\_       Puff  Plast  P.Cell  Aliplast
  - Calf        \_\_\_\_\_       Puff  Plast  P.Cell  Aliplast
  - Other: \_\_\_\_\_ \_\_\_\_\_     Puff  Plast  P.Cell  Aliplast

- Quality Assurance**
- \_\_\_\_ Smooth edges      \_\_\_\_ Joint Alignment      \_\_\_\_ Mods
  - \_\_\_\_ Symmetrical trim    \_\_\_\_ Comp. Loctited      \_\_\_\_ Pull
  - \_\_\_\_ Pads glued/Incl.    \_\_\_\_ All parts Incl.      \_\_\_\_ Finish
  - \_\_\_\_ Straps On/Incl.      \_\_\_\_ Proj. tote cleared    \_\_\_\_ QA
  - \_\_\_\_ Strap placement    \_\_\_\_ Engraved P.O.      \_\_\_\_ Finish Date

- Modifications**
- Standard modifications are no reductions on AFOs with 1/8" ankle buildup

- Strut Selection**      Standard 10.5" x 1.5" inches / Mini 8.5"x1.25"
- |   |   |
|---|---|
| <b>Standard Dynamic Strut</b>                             | <b>Mini Dynamic Strut</b>                                 |
| <input type="checkbox"/> CD207XS - Extra Soft Flex Strut  | <input type="checkbox"/> CD209S - Soft Flex Strut         |
| <input type="checkbox"/> CD207S - Soft Flex Strut         | <input type="checkbox"/> CD209M - Medium Flex Strut       |
| <input type="checkbox"/> CD207M - Medium Flex Strut       | <input type="checkbox"/> CD209R - Rigid Flex Strut        |
| <input type="checkbox"/> CD207R - Rigid Flex Strut        | <input type="checkbox"/> CD209XR - Extra Rigid Flex Strut |
| <input type="checkbox"/> CD207XR - Extra Rigid Flex Strut |   |

- Anterior Shell**
- No Shell (Strap Only)
  - FASS - Simple Shell
  - FAAS - Molded Shell

- Insole Support / External Posting**
- ZPIS Pelite Insole (per layer)
  - Neutral to Bench
  - Distance from Neutral \_\_\_\_\_

**Hourly Labor**     FHL \_\_\_\_\_

**Mod Notes**

---

---

---

---

---

---

---

---

**Additional Notes**

---

---

---

---

---

---

---

---

Send cast/parts to: Coyote, 525 N. Steelhead Boise, ID 83704